


1MEDiX Provider Portal Screen Mockups

Navigation Menu

Navigation Menu expanded to show the following sections:


1. Dashboard
2. Calendar - Consolidated calendar view
3. Messages - Secure message / DIRECT message center
4. Contacts - Displays contacts and providers in the current organization
5. Patients - List of the patients associated with this provider. If patients opt out from the provider, they will not be on this list.
6. Reports - Query and reporting tools
7. Admin - Administrative tools



1MEDiX Provider Portal Screen Mockups

Dashboard

1. Navigation menu
2. Message Center alerts - displays number of new messages
3. Displays user alerts
4. User menu - user profile, user preferences, user security settings
5. Customizable and moveable Dashboard widgets - widgets are customized per deployment needs



The dashboard mockup features a blue header bar with the 1MEDiX logo and navigation menu (1). On the right, there are alert icons: a message icon with '5' (2), a bell icon with '0' (3), and a user profile dropdown for 'Dr. Douglas Howser, M.D.' (4). Below the header, the dashboard is divided into two main sections. The left section, titled 'New Messages' (5), lists five messages from 'United Hospitals of America' with dates ranging from 1/1/16 to 1/11/16. Each message includes a 'LAB' icon and a title. The right section, titled 'Patients', displays a table of patient information. Below the table is a 'View More Patients' button. The 'New Messages' section also includes a 'View More Messages' button at the bottom.


Patient	Gender	Date of Birth
Judy Glenn	female	Mar 30, 1986 (30 yr. 1 mo.)
Tiffany Westin	female	May 20, 1975 (41 yr.)
Pedro Burton	male	Jul 4, 1952 (63 yr. 10 mo.)
Bruce Little	male	Jan 4, 1943 (73 yr. 4 mo.)
Dustin Gomez	male	Mar 15, 1949 (67 yr. 2 mo.)

Message Center

1. Category - icons indicate the message category: Laboratory, Radiology, etc.

1MEDiX Provider Portal Screen Mockups



2. From - Message source (person or facility)
3. Subject - Current example is a message created from a HL7v2 ORU^R01
4. Patient - patient associated with the message
5. Date - message sent date





1MEDiX Provider Portal Screen Mockups

Message Content

1. Message information
2. Star - mark the message as a favorite for easy finding
3. FHIR Resources attached to message - Current example is Diagnostic Report from HL7v2 ORU^R01 Lab report



4

0

Dr. Douglas Howser, M.D. ▾

Messages > CBC Differential Auto, CBC Hemogram

CBC Differential Auto, CBC Hemogram

Laboratory Result

2


☆ Star

1

From United Hospitals of America
Patient [Tiffany Westin](#)
Date 1/11/16 8:06 AM

3

CBC Differential Auto
Effective Date 1/11/16 7:26 AM
Issued Date 1/11/16 8:06 AM
Status final

Result	Value		Reference Range	Comments
Neut%: Neut%	68.4 %		40.0-74.0 (%)	
Lymph%: Lymph%	23.4 %		19.0-48.0 (%)	
Mono%: Mono%	6.1 %		3.5-9.0 (%)	
Eos%: Eos%	1.5 %		0.0-7.0 (%)	
Baso%: Baso%	0.7 %		0.0-1.5 (%)	

1MEDiX Provider Portal Screen Mockups

Patient Profile - Basic Info

1. Header - basic information about patient
2. Export - export the current patient into various formats (CDA, FHIR, etc)
3. Contact - contact the patient by secure message, phone, etc
4. Profile Tabs - customizable tabular display of patient data
5. Basic demographic information of patient

1MEDiX

Patients > Tiffany Westin > Basic Info

Tiffany Westin
DOB: May 20, 1975 | 41 yr.
Gender: female

Export Contact Patient

Basic Info Allergies Care Team Conditions Diagnostic Reports Documents
Encounters Immunizations Medications Observations Procedures

Identifiers Medical record number 8888802

Names official Tiffany Westin

Birth Date May 20, 1975 (41 yr.)

Contact Info Personal (626) 555-6734


Addresses Personal 128 Dempsey Rd
Columbus OH 43213


Patient Profile - Diagnostic Reports


1. Diagnostic Report information
2. Button to show/hide associated Results (Observation)

1MEDiX Provider Portal Screen Mockups

3. Result (Observation) information
4. Button to show/hide Observation charting
5. Each column is sortable

 1MEDiX

 4

 0

Dr. Douglas Howser, M.D. ▾

Patients > Tiffany Westin > Diagnostic Reports

**Tiffany Westin**DOB: May 20, 1975 | 41 yr.
Gender: female

Export ▾

Contact Patient ▾

Basic Info

Allergies

Care Team

Conditions

Diagnostic Reports

Documents

Encounters

Immunizations

Medications

Observations

Procedures


Show 10 entries

Search:

1MEDiX Provider Portal Screen Mockups

Patient Profile - Diagnostic Reports - Chart Expanded

- 1. High range is shaded in red
- 2. Normal range is shaded in green
- 3. Historical value points (clicking/hovering over points will display more information)



1MEDiX Provider Portal Screen Mockups


Patient Profile - Mobile view

1. Patient Tabs condensed into a dropdown menu
2. Responsive table

1MEDiX

Dr. Douglas Howser, M.D. ▾

Patients > Tiffany Westin > Conditions

 **Tiffany Westin**
DOB: May 20, 1975 | 41 yr.
Gender: female

1 ▾ Conditions ▾ Contact Patient ▾ Export ▾

Show 10 entries

Search:

Date	Condition	Status
1/1/85	195967001: Asthma	confirmed
1/1/85	61582004: Allergic rhinitis	confirmed
1/1/95	89765005: Tobacco use disorder	confirmed

1MEDiX Provider Portal Screen Mockups

Patient Profile - Mobile view - Expanded

The Patient Profile tabs are expanded showing a drop-down menu

